6 MeV X-Band Li

6 MeV X-Band Linac

NDT & Medical Applications

J. Mondal, L. Mishra, S. G. Sarkar, V. Sharma, V. Kaushik, R. Kalra, N. Chaudhary, E. Kandaswamy, R. I. Bakhtsingh and Archana Sharma

band low energy accelerators are compact in size, light weight, which make them highly superior over S and C band Linacs, with applications in medical radiation therapy/non destructive testing. BARC has developed a 9300 MHz, 6 MeV, 0.48 kW standing wave type X-Band Linac (XBL) cavity based x-ray source. The X-Band Linac based x-ray source has been assembled with primary collimator and dose measurement setup at a distance of 1 m from the source. The 6 MeV XBL comprises an on-axis coupled cavity configuration with 49 cells to it(25 accelerating and 24 coupling cells). The Linac is operated in the pulsed mode with a pulse width of 3-4 μs at 200 Hz repetition rate. An indirectly heated thermionic dispenser cathode based electron gun produces pulsed electron beam, which is then injected into the Linac cavity at 20-25 keV energy. The XBL is fed with RF power from a magnetron based power source at 9300±10 MHz to generate an accelerating electric field gradient of approximately 18.5 MV/m in vacuum of 1×10^{-6} mbar or higher. Post acceleration, the 6 MeV electron beam is focused to less than 2 mm diameter with the help of a permanent focusing magnet. The focused beam falls on a water-cooled Tungsten target to produce x-rays of 600 cGy/min at a distance of 0.8 m from the target at 200 Hz. This x-ray beam is then collimated with a 12.7° conical collimator to obtain a circular x-ray beam of diameter of

The authors are from Accelerator & Pulse Power Division,

Laser & Plasma Technology Division and Electron Beam Centre of Beam Technology Development Group, BARC.

around 97 mm at the exit of the collimator. PLC-based control and interlock (Control & Instrumentation) system of the Linac based x-ray source has been implemented for fail-safe operation. In the event of any interlock failure, both the primary electron beam and RF source power are switched OFF to maintain the integrity of the system.

The authors gratefully acknowledge the active role and support provided by CDM, BARC and EBC, Kharghar during the development of Linac.

Applications

Inter-Operative Radiation Therapy

Stereotactic Radiosurgery

Non Destructive Testing of materials and industrial components