

Simultaneous Solid State Detector ICP-OES (Activa)

FEATURES

- ❑ Wavelength range from 165-800 nm
- ❑ Grating of 4343 and 2400 grooves per mm
- ❑ Charge Coupled Device (CCD) solid state detector.
- ❑ Solution uptake 0.9 ml/min


ADVANTAGES

- ❑ 75 elements
- ❑ Resolution up to 10 picometer
- ❑ 6 mm height of the plasma view
- ❑ Signal to Noise Normal Analytical Zone view
- ❑ Fast



Sample Analysis

- Ores, minerals
- Alloys
- Solutions



[back](#)